


# BRIEF INSTRUCTIONS

## AT-F7-HV



**JJ X-RAY**  
Danish Science Design



## AT-F7-HV


This slit is a further development of the original ESRF slit (AT-F7-AIR) which was only designed for application in air or inside a larger vacuum chamber. For this slit system, we have extended the ESRF design to allow for connections via vacuum flanges.

The characteristics of this system is the same as the regular ESRF slit: The compact size, the sub-micron resolution and the blade control (where one motor defines the aperture size and the other one determines the position of the aperture).



Please note the aperture size and the relatively low motor current.

# Technical Specifications

AT-F7-HV	
Aperture size	Maximum: 7 mm x 7 mm Minimum: Full overlap
Resolution	Aperture: 0.813 micron per full step Translation: 0.322 micron per full step
Repeatability	< 1 micron
Accuracy	± 3 micron (over 3 mm)
Vacuum	O-ring sealed, high vacuum 10-5 mbar, low outgassing materials
Mechanical dimension	145 mm x 145 mm x 80 mm
Standard blades	2 mm thick tungsten carbide with knife-edge profile (2 degrees slope)
Mechanical connections	M6-threaded holes on the sides as shown on the drawing. KF-40 flange connections are provided on the body of both sides of the slit
Limit switches (end-of-travel)	Included as standard on all motions
Weight	≈ 1.6 kg
Outer surface	Anodized aluminum in color nature
Guiding	High precision internal rails and carriages
Electrical connections	Microswitches coupled to 2 motor connectors (SUB-D 15 pins male)
Motors	2 phase stepping motors



Issue	29-01-2003	File		Part	
Client	JJ X-RAY AS	Project	JJ X-Ray AT7HV Silt	Scale	A3
With line source		Name		Complete	AT7HV



JJ X.RAY Michael Andersen	
File: A17airrv11.sch	
Sheet: /	
Title: A1-F7-AIR/HV	
Size: A4	Date: 19 Jan 2018
KiCad E.D.A.	eesschema (2013-05-16 BZR 4016)-stable
Rev 1.1	Id: 1/1

# Motor Specifications

Motor Specifications	
Number of motors	4
Motor type	2-phase stepping motor
Step angle	15°
Connection type	Bipolar (Serial)
Current per phase	0.45 A/phase
Resistance	3.6 Ω/phase
Inductance	1.9 mH/phase
Limit switches	'+' and '-' end of travel

Motion Mechanism, Aperture (M1)	
Type of motion	Two translations in opposite directions
Guidance	In vacuum linear rails and carriages
Motor step angle	15°/step
Motor gear	41:1
Lead screw pitch	2x0.4 mm/rev (left and right)
Scale factor, aperture	1230 steps/mm (gap)
Mechanical resolution	≈0.813 μm/step
Aperture calibration	≈1230 steps/mm

Motion Mechanism, Translation (M2)	
Type of motion	Translation
Guidance	In vacuum linear rails and carriages
Motor step angle	15°/step
Motor gear	41:1
Pitch	0.3175 mm/rev
Scale factor, translation	3099.2 steps/mm
Mechanical resolution	≈0.3227 μm/step
Translation calibration	≈3099.2 steps/mm

Recommended Driver Settings	
The motors should be run <u>at 0.45 A per phase.</u> The motors have been tested at:	
Running speeds	800-1000 steps/second
Starting speeds	300-400 steps/second
Ramp times	0,1 second

Always use “backlash correction” if available (i.e. the motor always approaches the final position from the same side). A useful backlash parameter could be 0.1 mm.

Warning:


If you are using systems/detectors that can be damaged by overexposure, where the slits are used to remove a lot of the intensity, be careful when changing aperture size since the backlash correction may result in the slit being opened significantly more than you anticipated during adjustment.

## Manual Control

If for some reason, you need to move the slit blades manually, it is possible to use the screw attached to the back shaft of the motor. This can be done by inserting an M3 Allen wrench into the slit housing to engage the screw. It is probably easiest if this operation is done with an open cover so you can see what is going on (firstly: Switch off the controller). In the vacuum slit versions, you will need to rotate in increments using the Allen wrench.

## Limit Switches

The limit switches should be wired up, if at all possible. There is always some ambiguity in the definition of the travel direction. The first time you test the slits and cables you should therefore open up the slit and test the actual functioning of the limit switches. Below we show an image that may help you in determining the appropriate limit switch setup.



*Mechanical overview of limit switch position  
M1 movement: Aperture  
M2 movement: Translation*

# Trouble Shooting

The most common issues and their resolutions are:

- ❏ The motor does not move when it should (it can be silent, be jittering or be making a noise).
  - The motor is not receiving enough current. Try setting the current a bit higher (for example 10%). If problems persist check with an amp-meter to see that your driver is working properly.
  - The wiring is bad. Check cabling.
  - One of the motors' phases is burnt. Check that the resistance on all phases is the same. If not, contact us to have the slit sent for repair.
  - The motor is stuck against a limit switch. Un-stick it, using the scale wheel or, if an AT-slit, open the slit (see manual control), and fix the limit switch issue.
- ❏ Restart the controller and the controller program.
- ❏ The blade system shows irreproducibility during operation.
  - The rail system may have become loose. Open the slit. Check if the rail-system is tight. Tighten screws if you need to

## Common options

- ❏ Blade options: 4, 5 or 10 mm, other blade materials.
- ❏ Connectors: Other KF- and CF-flanges.
- ❏ Special preparation for enabling vacuum  $<10^{-6}$  mbar.



# The JJ X-Ray Product Range

- ☒ Slit Systems (AIR, HV, UHV)
- ☒ Complete Beamline Solutions
- ☒ Spectrometers
- ☒ Refractive Optics
- ☒ Foil Collimators
- ☒ Positioning

## Contact JJ X-Ray A/S

If you have any questions, concerns, request for quotations or need general advice, please feel free to contact us:

**JJ X-Ray A/S**  
Scion-DTU, Dr. Neergaards Vej 5D  
2970 Hoersholm  
Denmark

Telephone: +45 4776 3000  
Fax: +45 4776 3001  
VAT: DK 29523215

info@jjxray.dk  
sales@jjxray.dk  
Website: [www.jjxray.dk](http://www.jjxray.dk)